

SHARPSGUARD® *eco yellow 2.5*

Product Code: DD472YLR

Our SHARPSGUARD® *eco yellow* range is used for the disposal of sharps that have been contaminated with medicines.

This container is manufactured using 44% recycled material and is certified to ISO23907:2019 and UN3291.

This container has a yellow base and a range of coloured lids so is compliant with HTM 07-01 guidance for the safe management of healthcare waste.

Our SHARPSGUARD® *eco yellow 2.5* is a general-purpose sharps container, it provides the maximum capacity in the minimum working space.

- Syringe Sizes: 1, 2, 5, 10, 20ml
- No. of Syringes*: Up to 120
- Assembly: Lid and base supplied separately
- Accessories available:
 - DD487 POUDS® small wall bracket
 - DD497 POUDS® round leg trolley bracket
 - DD501 POUDS® square leg trolley bracket
 - DD505 POUDS® medi rail bracket
 - PT199D POUDS® combination tray
 - PT299 POUDS® tray 2.5

*based on 1ml Insulin needle with a 29g needle attached



Usable Capacity (l)	Product Weight (g)	Product Dimensions Assembled (max excluding bracket handles)			Aperture Dimensions (max)		Handle
		Height (mm)	Width (mm)	Diameter / Depth (mm)	Width (mm)	Depth (mm)	
2.21	232	166	N/A	176	77	46	Yes

Carton Quantity	Carton Dimensions					Carton Weight & Containers (Kgs)
	Height (mm)	Width (mm)	Depth (mm)	Carton Weight Empty (g)		
48	346	663	343	550		11.68



Safety, Service and Support – *Every Day, Affordably & Responsibly*

www.daniels.co.uk

